Form 100 - Blood Collection and Processing

OMB # 0925-0414 Exp: 4/06

	- Affix label here-	
Clinical Center/II	D:	
First Name		M.I
Last Name		<u> </u>

BLOOD REQUEST (Mark one.) ____HRT ___DM ___HRT + DM ___OS

Blood		Three 7 ml		Two 4.5 ml		One 10 ml		One 2 ml	
Collection Tube		Royal Blue		Light Blue		Lavender			Lavender
			Trig.	Coag Panel	One 1.8	Lipid Panel			
		Four 1.8 ml	0.5 ml	Three 1.8	Buffy	Three 1.8		1 Buffy	
Ст	ryovials	Serum	Serum	ml Plasma	coat	ml Plasma	1 RBC	Coat	CBC
Study	Visit								
CT/OS	SV1	Х	HRT if Lipemic	Х	Х	Х	Х	Х	Х
СТ	1st Annual	х	HRT if Lipemic	х	Х	Х	Х	Х	
СТ	Subsample at 3rd, 6th, and 9th Annual	Х	HRT if Lipemic	Х		Х			
OS	3 Year	Х		Х	Х	Х	Х	Х	Х

BLOOD COLLECTION

1.	Date blood drawn:		(M/D/Y)				
2.	Drawn by:						
3.	Contact type:	□_ ₁ Phone	□_ ₃ Visit				
		Mail	☐ ₈ Other				
4.	Visit type:	□_1 Screening	#				
		Semi-Annual	[#] ∟I				
			#				
		\square_4 Non-Routine					
5.	Time drawn:	(Hr:Min)	□ ₁ AM □ ₂ PM				
BLOOI	O COLLECTION CHECKLIS	т					
	"How long since you had anything to eat or drink besides water?" hours (If you are drawing for a lipid panel and this is less than 12 hours, do not draw blood. Arrange a time when the woman can come in for a fasting blood draw.)						
6.	(If you are drawing for a lipi	d panel and this is less than 12	hours, do not draw	hours			
6. 7.1.	(If you are drawing for a lipi blood. Arrange a time when	d panel and this is less than 12	hours, do not draw asting blood draw.)	\square_0 No \square_1 Yes			
	(If you are drawing for a lipit blood. Arrange a time when "Have you engaged in any	d panel and this is less than 12 n the woman can come in for a f	hours, do not draw asting blood draw.) the last eight hours?"				

К_____

BLOC	DD PROC	ESSING					
9.	Processed by:						
10.	Time began centrifugation:		•	(Hr:Min)	□ ₁ AM □ ₂ PM		
11.	For HRT and HRT + DM Screening Visits and Annual Visits Only: Is serum in royal blue tube lipemic? \Box_0 No \Box_1 Yes \longrightarrow Process aliquot for triglyceride level						
12.	Time s	ample placed in cryovi	als:	(Hr:Min)	□ ₁ AM □ ₂ PM		
13.	Time c	ryovials placed in freez	zer:	(Hr:Min)	□ ₁ AM □ ₂ PM		
14.	WHI bl	ood sample number:	- Affix blood sample "Form" label here -				
		Orig Tube Color	Cryovial Color/Test	15. Cryovial Number	16. Mark if Sample Processed		
		Royal Blue	Orange	0 2			
		Royal Blue	Orange	0 3	\square_1		
		Royal Blue	Orange	0 4	\square_1		
		Royal Blue	Orange	0 5	\square_1		
		Light Blue	Blue	<u>0</u> 6			
		Light Blue	Blue	0 7			
		Light Blue	Blue	0 8			
		Royal Blue 7ml	Trig (0.5 ml)	0 9			
		Lavender 10ml	Yellow	$1 \mid 0$	\square_1		
		Lavender 10ml	Yellow	1 1	\square_1		
		Lavender 10ml	Yellow	<u>1</u> 2	\square_1		
		Lavender 10ml	White (Buffy Coat)	1 3	\square_1		
		Lavender 10ml	Red (RBC)	<u>1</u> 4			
		Lavender 2 ml	CBC	<u>1</u> 6			
		Light Blue 10 ml	White (Buffy Coat)	2,0	\square_1		